SALTER

Body Analyser Scale

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BODY FAT – WHAT DOES IT MEAN?

The human body is made up of, amongst other things, a percentage of fat. This is vital for a healthy, functioning body – it cushions joints and protects vital organs, helps regulate body temperature, stores vitamins and helps the body sustain itself when food is scarce. However, too much body fat or indeed too little body fat can be damaging to your health. It is difficult to gauge how much body fat we have in our bodies simply by looking at ourselves in the mirror. Hence the need to be able to measure and monitor body fat percentage.

Body fat percentage gives you a better measure of fitness than weight alone – the composition of your weight loss could mean you are losing muscle mass rather than fat – you could still be ‘over-fat’ even when a scale indicates ‘normal weight’.

BODY WATER MEASUREMENT – WHY MEASURE IT?

Body water is the single most important component of body weight. It represents over half of your total weight and almost two thirds of your lean body mass (predominantly muscle).

Water performs a number of important roles in the body:

All the cells in the body, whether in the skin, glands, muscles, brain or anywhere else, can only function properly if they have enough water. Water also plays a vital part in regulating the body’s temperature balance, particularly through perspiration.

The combination of your weight and fat measurement could appear to be ‘normal’ but your fluid intake could be insufficient for healthy living.

HOW CAN KNOWING YOUR BODY FAT AND WATER PERCENTAGES HELP?

It can serve as a guide to whether you need to adjust your diet, fitness program or fluid intake to help achieve a healthy balance. Refer to tables (See Fig 1 and Fig 2).

HOW DOES THE BODY FAT ANALYSER SCALE WORK?

Fat cells and lean/muscle cells have different conductivity. Using a system known as Bio-electrical Impedance Analysis (BIA) the body fat analyser passes a minute, harmless, electrical impulse through the body to determine a resistance reading. This reading is then processed along with the users age, weight, sex and height to calculate the body fat and body water percentage. For the most consistent results only use the Salter Body Analyser Scale when your body is normally hydrated (i.e. not following excessive perspiration, nor just after a bath, nor after drinking a lot of water or alcohol).

QUESTIONS AND ANSWERS

How exactly is my body fat and water being measured?

The Salter Body Analyser Scale uses a measurement method known as Bioelectrical Impedance Analysis (BIA). A minute current is sent through your body, via your feet and legs. This current flows easily through the lean muscular tissue, which has a high fluid content, but not easily through fat. Therefore, by measuring your body's impedance (i.e. its resistance to the current), the quantity of muscle can be determined. From this, the quantity of fat and water can then be estimated.

What is the value of the current passing through me when the measurement is taken?

Is it safe?

The current is less than 1mA, which is tiny and perfectly safe. You will not be able to feel it. Please note however, that this device should not be used by anyone with an internal electronic medical device, such as a pacemaker, as a precaution against disruption to that device.

If I measure my body fat and water at different times during the day, it can vary quite considerably. Which value is correct?

Your body fat percentage reading varies with body water content, and the latter changes throughout the course of the day. There is no right or wrong time of the day to take a reading, but aim to take measurements at a regular time when you consider your body to be normally hydrated. Avoid taking readings after having a bath or sauna, following intensive exercise, or within 1-2 hours of drinking in quantity or eating a meal.

My friend has a Body Fat Analyser made by another manufacturer. When I used it I found that I got a different body fat reading to that on my Salter Body Analyser Scale. Why is this?

Different Body Fat Analysers take measurements around different parts of the body and use different mathematic algorithms to calculate the percentage of body fat. The best advice is to not make comparisons from one device to another, but to use the same device each time to monitor any change.

How do I interpret my body fat and water percentage readings?

Please refer to the Body Fat and Water tables included with the product. They will guide you as to whether your body fat and water reading falls into a healthy category (relative to your age and sex).

What should I do if my body fat reading is very 'high'?

A sensible diet, fluid intake and exercise program can reduce your body fat percentage. Professional medical guidance should always be sought before embarking on such a program.

Why are the body fat percentage ranges for men and women so different?

Women naturally carry a higher percentage of fat than men, because the make-up of the body is different being geared towards pregnancy, breastfeeding etc.

What should I do if my water percentage reading is 'low'?

Ensure that you are regularly taking sufficient fluids and work towards moving your fat percentage into the healthy range.

Why should I avoid using the Salter Body Analyser Scale whilst pregnant?

During pregnancy a woman's body composition changes considerably in order to support the developing child. Under these circumstances, body fat percentage readings could be inaccurate and misleading.

NORMAL HEALTHY RANGE OF BODY FAT & TOTAL BODY WATER PERCENTAGE

MALE			FEMALE		
AGE	BODY FAT	BODY WATER	AGE	BODY FAT	BODY WATER
6	12-20%	60%+	6	14-23%	50%+
7	13-20%		7	15-25%	
8	13-21%		8	15-26%	
9	13-22%		9	16-27%	
10	13-23%		10	16-28%	
11	13-23%		11	16-29%	
12	12-23%		12	16-29%	
13	12-22%		13	16-29%	
14	11-21%		14	16-30%	
15	10-21%		15	16-30%	
16	10-20%	60%	16	16-30%	50%
17	10-20%		17	15-30%	
18-29	8-19%	60%	18-29	20-30%	50%
30-39	11-21%		30-39	22-32%	
40-49	13-23%	55%	40-49	24-34%	47%
50-59	15-25%		50-59	26-36%	
60-69	17-27%	50%	60-69	28-38%	45%

Fig 1

	BELOW NORMAL BODY WATER %	NORMAL & ABOVE BODY WATER %
LOW BODY FAT %	Fitness program recommended, ensure sufficient water intake	Fitness program recommended
HEALTHY BODY FAT %	Weight and body fat OK but insufficient water intake	Ideal balance body fat, muscle mass and weight
HIGH BODY FAT %	Change of diet required, ensure sufficient water intake	Change of diet required

Fig 2

For those using Athlete Mode: note that athletes may have a lower body fat range and a higher body water range than shown above, depending on their type of sport or activity.

Body water measurement results are influenced by the proportion of body fat and muscle. If the proportion of body fat is high, or the proportion of muscle is low then the body water results will tend to be low.

It is important to remember that measurements such as body weight, body fat and body water are tools for you to use as part of your healthy lifestyle. As short term fluctuations are normal, we suggest you chart your progress over time, rather than focus on a single days reading.


This product is for domestic use only. Pregnant women should only use the weight function.

Professional medical guidance should always be sought before embarking on diet and exercise programmes.


CAUTION: Do not use this product if you have a pacemaker or other medical device fitted. If in doubt, consult your doctor.

PERSONAL DATA ENTRY


- 1) Press the **ON/OFF** button to turn the scale on.
- 2) **'HI'** will be displayed followed by **'1-12'**.
- 3) While **'1-12'** is flashing, press the **SET** button to activate the set-up procedure.
The current user number will flash.
- 4) While the user number is flashing, select a user number by pressing the **UP** or **DOWN** buttons.
Press the **SET** button to confirm your selection.

- 5) The male or female symbol will flash.
Press **UP** button until the correct symbol is flashing, then press the **SET** button.


- 6) The height display will flash.
Press **UP** or **DOWN** as necessary to set your height, then press the **SET** button.


- 7) The age display will flash.
Press **UP** and **DOWN** as necessary to set your age, then press the **SET** button.

Child Mode:
When the age is set from 6-12 years child mode is automatically selected.
A small child icon will appear, as illustrated below.

- 8) '1-12' will flash, then your height and gender will be displayed.
Wait for '0.0' to be displayed, then allow the scale to switch off. The memory is then set.

- 9) Repeat procedure for a second user, or to change user details.


Note: To update or overwrite the memorised data, follow the same procedure, making changes as required.

Tip: During the set-up procedure press **on/off** to quit and start again at any time.

WEIGHT READING ONLY

Position the scale on a flat level surface.

- 1) Press the **ON/OFF** button to turn the scale on.
- 2) 'HI' will be displayed followed by '1-12'.
- 3) Whilst '1-12' is flashing, press **ON/OFF** button to enter weight reading only mode.
- 4) Zero is displayed.
- 5) Whilst zero is displayed, step onto the platform and stand still.
- 6) After a few seconds your weight is displayed.

WEIGHT & BODY FAT READINGS

Position the scale on a flat level surface.

- 1) Press the **ON/OFF** button to turn the scale on.
- 2) 'HI' will be displayed followed by '1-12'.
- 3) While '1-12' is flashing, select your user number by pressing the **UP** or **DOWN** buttons.
Your height and gender will be displayed.
- 4) Zero will be displayed.


Athlete Mode:

To select Athlete Mode press the **UP** button. A small athlete icon will appear.

Note: This feature is only available to those aged between 15-50 years.

An athlete is defined as a person who is involved in intense physical activity of approximately 12 hours per week and who has a resting heart rate of approximately 60 beats per minute or less.

- 5) When zero is displayed, with bare feet, step onto the platform and stand still.
- 6) After 2-3 seconds your weight will be displayed.
- 7) Remain standing on the scale while your other readings are taken.
- 8) Your body fat % (F) and total body water % (W) are displayed.
The results are repeated then the scale switches off.

WARNING INDICATORS

'Lo' = replace battery.

'O-Ld' = Weight exceeds maximum capacity.


'Err' = out of range, incorrect operation or poor feet contact.

TECHNICAL SPECIFICATION

80 cm-220 cm
6-100 years
Male/Female
12 user memory


Child mode = 6-12 years
Athlete mode = 15-50 years
d = 0.1 % body fat
d = 0.1 % body water

BATTERY


Requires 2 x CR2032 batteries.

POSITION

Position scale on flat surface.

WEIGHT MODE

Set the switch on underside of product to: KG, ST, or LB.

ADVICE FOR USE AND CARE

- Always weigh yourself on the same scale placed on the same floor surface. Do not compare weight readings from one scale to another as some differences will exist due to manufacturing tolerances.
- Placing your scale on a hard, even floor will ensure the greatest accuracy and consistency.
- It is recommended that you measure at the same time of the day, preferably early evening before a meal for the most consistent results.
- Your scale rounds up or down to the nearest increment. If you weigh yourself twice and get two different readings, your weight lies between the two.
- Clean your scale with a damp cloth. Do not use chemical cleaning agents.
- Do not allow your scale to become saturated with water as this can damage the electronics.
- Treat your scale with care – it is a precision instrument. Do not drop it or jump on it.

TROUBLESHOOTING GUIDE

- You must have bare feet when using body fat and total body water measurements.
- The condition of your skin on the bottom of your feet can affect the reading. To get the most accurate and consistent reading, wipe your feet with a damp cloth, leaving them slightly moist before stepping on the scale.

WEEE EXPLANATION

 This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

GUARANTEE

This product is intended for domestic use only. Salter will repair or replace the product, or any part of this product, (excluding batteries) free of charge if within 15 years of the date of purchase, it can be shown to have failed through defective workmanship or materials. This guarantee covers working parts that affect the function of the scale. It does not cover cosmetic deterioration caused by fair wear and tear or damage caused by accident or misuse. Opening or taking apart the scale or its components will void the guarantee. Claims under guarantee must be supported by proof of purchase and be returned carriage paid to Salter (or local Salter appointed agent if outside the UK). Care should be taken in packing the scale so that it is not damaged while in transit. This undertaking is in addition to a consumer's statutory rights and does not affect those rights in any way. For UK Sales and Service contact HoMedics Group Ltd, PO Box 460, Tonbridge, Kent, TN9 9EW, UK. Helpline Tel No: (01732) 360783. Outside the UK contact your local Salter appointed agent.

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